

Identification and Risk Mitigation of a Critical Process Parameter during Antibody Process Scale-up

Christa Short, Midge Bonzon, Jonathan Mitschelen, Upstream Process Development, R&D Services (Reference)

Background

Start-up biotechnology companies are founded on great product concepts from partners or universities

- Early partners and investors want to see proof of concept
- Funding is often in tranches tied to development milestones

Entrepreneurial Business Strategies

- Rapid realization of value
- Limited financing
- Limited duration commitment-early exit preferred
- Minimal business organization or contract experience
- High risk tolerance – “Plan for Success”
- Plans change, assumptions do not

Impact on Manufacturing

Timeline

- Development of manufacturing process is rate-limiting activity for progression to clinical studies
- Virtual or small company needs 18 – 24 months to develop process
- Investor timelines often require clinical trial initiation in 12-14 months

Process Development Strategies

- Process and Process Knowledge “Adequate” for Business Plan
- Graduated cGMP development and compliance
- Delay major development expenses - “Fix before Phase III”
- Rely on minimally characterized standard processes
- At risk scale up – Straight to cGMP Production for Phase II

Company Strategies

To meet timelines and reduce costs prior to proof of concept, companies often cut corners

- Processes deliver safe, effective products at small scale
- Process parameters and process design space are unmapped
- Understanding likely cell culture behavior during scale-up based on scale-down models is not a priority

Material is available for early clinical trials, but...

- Patient safety is never at risk, nor is mode of action
- Ability to meet ongoing clinical demand IS at risk due to limited process knowledge
- Rapid requirement to scale-up to meet clinical supply can lead to unexpected surprises...even with antibodies

Regulatory Consideration

FDA Phase 1 Guidance

- Phase appropriate GMPs
- Emphasis on patient safety but limited expectation for process consistency early in development
- Continuous improvement in process control throughout development
- Operation within historical ranges not considered a process change

GMPs for the 21st Century

- Continuous process improvement
- Knowledge driving process design
- Risk-based decision making driven by data and scientific rationale

Case Study

- Production of Phase 1 clinical trial material at 100 L scale
 - Limited process development prior to cGMP production
 - Productivity of < 1 g/L, below industry standard
- Scale-up to 1,000 L planned for Phase 2 clinical supply
 - Initial 1,000 L batch successful with similar yield to 100 L batches
 - Second 1,000 L batch failed to reach normal cell density and yield was much lower than normal
- Risk mitigation for third batch included rigorous data review and implementation of a process control strategy

Mab Process Flowchart





Figure 1. PD & MFG Growth & Production.




Figure 3. Process Parameter Correlations.

Glutamate Issues & Process Range

Justification of Glutamate as a Process Issue

- Glu is a component of both base medium and feed medium
- Glu vs Titer strongest of correlations
- Glu has roles in both protein synthesis and energy production
- Literature - Metabolomics report Asn, Asp, Glu, Pyr depletion at start of decline phase associated with interruption of TCA cycle
- Schedule did not allow time for AA depletion studies
- Schedule did not allow time to establish Asn, Asp, Pyr assays

Process Range Considerations

- Consistent with regulatory considerations
- Lower action level reliably distinguishable from zero – 1.5 mM
- Range wider than maximum daily consumption – 2 mM / day
- Upper level: 1.5 mM + 2 mM ≈ 3.5 mM vs historical max = 3 mM

Mab Production Process

Product: Fully human IgG1
 Cell Line: CHO System
 Growth Medium: Gibco OptiCHO™
 Production Medium: Gibco OptiCHO™
 10% Gibco EfficientFeed™ A
 10% Gibco EfficientFeed™ B
 50% EfficientFeed™ A / 50% EfficientFeed™ B
 5% v/v of Initial Culture Vol/day Day 2 to Harvest
 Glucose as required
 Temperature: 37°C
 pH: 7.05 ± 0.15
 Harvest: ≤ 40% Cell Viability or Day 14

Process Performance History

- Process Development**
 - Problematic
 - Variable peak VCD
 - Variable run duration
 - Modest productivity
 - High oxygen consumption
- Manufacturing**
 - All lots met product specifications
 - 100 L & 200 L Scale
 - 5 runs
 - Yield ≥ PD runs
 - Variable duration, VCD, and yield
 - 1000 L Scale Run 1
 - Low vcd
 - Good run duration
 - Low but adequate yield
 - 1000 L Scale Run 2
 - Shortest run duration
 - Lowest, inadequate yield




Figure 2. PD & MFG Metabolite Profiles.

Process Change Options

- No Change**
 - Justify as maintenance of lot to lot process consistency
 - Inadequate data to support process change impact
- Supplement with Glutamate**
 - Justify as continuous process improvement
 - As required to maintain a standard control window
 - As required to maintain to range in successful runs
 - Adjust to standard concentration once at bioreactor inoculation
 - Adjust to standard concentration once at N-1 bioreactor inoculation

Final Change and Justification

Convert Process Options to Approved Manufacturing Change

- Agreement among Client, Regulatory & Process Consultants & CMO
- Reasonable manufacturing procedure
- Minimize potential for impact on product quality
- Minimal change to improve process control
 - Correct process deficiency: glutamate exhaustion
 - Target average Glu to historical high concentration: ~3 mM Glu

Agreed Process Changes

- Supplement N-1 Seed Train at inoculation only: 2.5 mM Glu
- Maintain Glu in production bioreactor: 2.5 – 4.5 mM Glu

Results




Figure 4. 1K L Process Change Culture Properties.




Figure 5. 1K L Process Change Metabolites.

- Lot met release specifications
- Yield satisfied critical clinical inventory requirements
- No process parameters significantly outside process history
- Most process parameters more favorable than historical range
 - Yield better than first 1K bioreactor run
 - Peak VCD in process history range
 - Duration to 14 days with viability >60%
- Changes implemented without introducing dramatic changes

Conclusions

- Improved Glutamate control successfully addressed scale up issue
- Intensive, collaborative process data 28.9
- review identified process change options that satisfied client business constraints, regulatory obligations, and manufacturing facility design

Reference

This poster is a version of slides presented at PEGS Summit in Boston, MA (May 5-9, 2014) in the session of “Scaling Up and Down Strategies for Protein Production”. The presentation was a joint effort with Susan Dana Jones, Ph.D., Vice President and Senior Consultant, BioProcess Technology Consultants.

Acknowledgements

Manufacturing, QA, QC, Project Management, BD & Client